## **HI-FI Temp**

Termómetro infrarrojo  
Infrared thermometer  
Thermomètre infrarouge  
Infrarot Thermometer  
Termometro a infrarossi  
Termômetro infravermelho  
Infrarood thermometer  
Termometr na podczerwień  
Υπέριθρο θερμόμετρο  
Инфракрасный термометр  
Termometru cu infraroșu  
Инфрaчeрвeн тeрмoмeтър  
ميزان حرارة

# Español

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**Fábrica:** Zhongshan Jinli Electronic Weighing Equipment Co., Ltd.

283rd South Min'an Road, Xiaolan Town,  
Zhongshan City, Guangdong, China  
Tel.: +86(760)28133810/28133793  
+867602813377828133779

**Representante CE:** Wellkang Ltd

Suite B,29 Harley Street  
LONDRES, W1G 9QR,U.K.  
www.CE-Marking.EU

# English

## Infrared thermometer HI-FI Temp

### Dear customer,

Many thanks for choosing to purchase a Taurus brand product.

Thanks to its technology, design and operation and the fact that it exceeds the strictest quality standards, a fully satisfactory use and long product life can be assured.

Read these instructions carefully before switching on the appliance and keep them for future reference. Failure to follow and observe these instructions could lead to an accident.

### SAFETY ADVICE AND WARNINGS

Use of this thermometer is not intended as a substitute for consultation with your physician.

There may be a risk of exasperating the state of an illness by self-estimate.

Please do not put our products close to the charged object in order to avoid the possibility of electric shock.

The user should contact the manufacturer if require further advice for using or maintaining your device or to report unexpected operation or events.

It is necessary to clean the device between uses on different PATIENTS, including methods for rinsing, drying, handling and storage until re-use.

Don't expose this thermometer to extreme temperature conditions of  $> 50\text{ }^{\circ}\text{C}$  or  $< -20\text{ }^{\circ}\text{C}$ .

Do not apply a strong shock to, drop, step on, or vibrate the main unit.

Do not drop or knock the device, and do not use if damaged.

Don't use the device beyond intended use.

Keep the device away from water and heat, including direct sunlight.

Do not use a portable phone near the unit. ·

Do not disassemble, repair, or modify the unit. ·

Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit. ´

Meanings of Symbols	
Signs	Notes on the signs
	Type BF Applied Part
	Refer to instruction manual
	Manufacturer
	European community
	Specific batch code
	CE Mark
	This marking shown on the product or its literature, indicates that it should not be disposed of, with other household wastes at the end of its working life. To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate this from other types of wastes and recycle it responsibly to promote the sustainable reuse of material resources.
<b>IP22</b>	IP22: The first number 2: Protected against access to hazardous parts with a finger, and the jointed test finger of 12 mm $\Phi$ , 80 mm length, shall have adequate clearance from hazardous parts and protected against solid foreign objects of 12,5 mm $\Phi$ and greater. The second number: Protected against vertically falling water drops when enclosure tilted up to 15°. Vertically falling drops shall have no harmful effects

	when the enclosure is tilted at any angle up to 15° on either side of the vertical.
--	---

3-10 years	36.1 to 37.8	97.0 to 100.0
11-65 years	35.9 to 37.6	96.6 to 99.7
> 65 years	35.8 to 37.5	96.4 to 99.5

### Technical data

<b>Normal operating conditions</b>	Temperature: 5°C~40°C
	Relative humidity: 15%~95%RH
	Atmospheric Pressure :860-1060hPa
	Power supply: DC 3V (2 x AA LR6)

### Intended Use

The forehead thermometer is designed for body surface and forehead temperature measurement for infant and adults without contact to human body.

The forehead thermometer can also be used to measure the temperature of a baby-bottle or bath or room temperature (by using the surface mode).

### Normal Temperatures According to measurement Method:

Measurement Method	Normal Temp °C	Normal Temp °F
Rectal	36.6 to 38	97.8 to 100.4
Oral	35.5 to 37.5	95.9 to 99.5
Axillary	34.7 to 37.3	94.4 to 99.1
Ear	35.8 to 38	96.4 to 100.4

The temperature of the human body varies throughout the day. It can also be influenced by numerous external factors: age, sex, type and thickness of skin...

### Normal Temperatures According to age:

Age	Temp °C	Temp °F
0-2 years	36.4 to 38.0	97.5 to 100.4

### Operation and setting method

#### Installing or replacing the batteries

- A. Push the battery cover back, open the cover from the body.
- B. Insert two "AA" batteries as indicated in the battery compartment and then replace the battery cover.

#### Setting method

Keep pressing "SET" key and hold for 3 seconds into the setting mode.

#### Select measurement unit (Menu F1)

Enter into setting mode, the LCD displays "F1", press "-" key to select "F" unit, press "+" key to select "C" unit.

#### Warning value setting (Menu F2) Default 37.8°C (100.4°F)

Enter into setting mode, short press "SET" key twice, the LCD displays "F2", then press "+" key to increase the setting value 0.1°C (0.1°F), press "-" key to reduce the setting value 0.1°C (0.1°F).

#### Offset value setting (Menu F3) Default 0.0°C (0.0°F)

Enter into setting mode, press "SET" key triple, the LCD display "F3", then press "+" key to increase the setting value 0.1°C (0.1°F), press "-" key to reduce the setting value 0.1°C (0.1°F).

#### Sound notification setting (Menu F4)

Enter into setting mode, press "SET" key four times, the LCD display "F4", then press "+" key to turn on the sound notification function, press "-" key to turn off the sound notification function.

#### Backlight setting (Menu F5)

Enter into setting mode, press "SET" key five times, the LCD display "F5", then press "+" key to turn on the backlight function, press "-" key to turn off the backlight function.

#### Exit setting mode

Press "SET" key until shutdown.

### Advice for the 1st Use

For stable and reliable results, it is essential to check the forehead thermometer and change as needed, as follows:

- Take the temperature of the same person using a conventional thermometer; you will get 37.0°C (98.6°F) for instance.
- Take the temperature of the same person using the forehead thermometer keeping the 50mm distance between the thermometer and the

forehead (Take care to remove any obstacle which could alter the measurement (hair, perspiration...). If you get 37.0°C (98.6°F) , the Non-contact body infrared Thermometer is properly set and ready for use.

If you get a lower temperature, 36.4°C (97.4°F ) for example, your difference is 0.6°C (1.2°F) , you should adjust the temperature on the forehead Thermometer Body infrared Thermometer and add the difference, i.e. 0.6°C ( 1.2°F) .

To do it, press the "SET" button for 3 seconds, the screen displays F1, press "SET" button again until you get F3, press "+" button in order to add the difference (in our example, 0.6°C (1.2°F) .

9.3.3 To check, take the temperature again using the Non-contact body infrared Thermometer.

### **Caution (this setting applies only to body mode)**

#### **Correct measurement**

##### **Measurement body temperature**

Hold the forehead thermometer, leave one finger in front of the button, keep transducer window at human's forehead, leaving a distance of 30-50mm from the forehead. Press the button, and the thermometer will start up automatically and display the measurement result.

##### **Measurement ambient temperature.**

Hold the thermometer, leave one finger in front of the button, Level the detection window face to the place where you want to measure, and then press button. The thermometer will start up automatically and display the measurement result.

#### **Cautions:**

- After replaced the batteries or move the thermometer to a new place, it is suggested to lay the forehead thermometer aside at least 30 minutes before using, so that you can measure temperature more accurately.
- If the measurement is made without an effective distance or a deflection of the central position of the measured object, a deflection from actual temperature may cause. It is suggested to repeat the measurement once or more.
- While measuring, it is suggested to focus the detection window at the forehead, and press the measurement button in front of the handle, scan the forehead around, so that it can measure the highest and lowest temperature of the forehead.
- It may affect the accuracy of measurements when the forehead is covered by hair, perspiration, cap or scarf.

### **Selecting measuring body temperature or surface temperature**

Turn on the thermometer, press "SET" key to select body temperature measurement mode or surface temperature measurement mode.

#### **Data Memory**

In measuring mode, press "+" or "-" key to view measurement result. "+" key for page up, "-" key for page down. Press the "+" or "-" key and hold for 2 seconds, the result switches quickly.

**Factory:** Zhongshan Jinli Electronic Weighing Equipment Co., Ltd.

283rd South Min'an Road, Xiaolan Town, Zhongshan City, Guangdong, China

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## English

### **WARRANTY AND TECHNICAL ASSISTANCE**

This product enjoys the recognition and protection of the legal guarantee in accordance with current legislation. To enforce your rights or interests you must go to any of our official technical assistance services.

You can find the closest one by accessing the following web link: <http://taurus-home.com/>

You can also request related information by contacting us.

You can download this instruction manual and its updates at <http://taurus-home.com/>

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