taurus

INCEPTION PERFECT

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Bathroom Scale

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Badwaage

Bilancia Pesapersona

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Badkamerweegschaal

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Bath scale INCEPTION PERFECT

DESCRIPTION

A	Platform
B	LCD screen
C	Setting buttons
D	Up / Down buttons
E	Male symbol
F	Female symbol
G	Age symbol
H	Height symbol
I	Weight unit (Kg, Lb, St)
J	Meaning off at measurement
K	Water measurement symbol
L	Muscle measurement symbol
M	Bone measurement symbol
N	Basal metabolic index
O	Battery compartment cover
P	Weight unit selector button.
Q	Electrodes
R	Health level indicator
S	Body Mass Index (BMI)

USE AND CARE

- Do not use the appliance if the on/off switch does not work.
- Do not overload the appliance's weight capacity.
- Remove the batteries from the appliance if you are not going to use it for some time.
- This appliance is for household use only, not professional, industrial use.
- This appliance isn't valid for trade.
- Store this appliance out of reach of children and/or persons with physical, sensory or reduced mental or lack of experience and knowledge.
- Do not open the battery/ies under any circumstances whatsoever.
- When the battery/ies is not in use, keep it away from other metal objects like paper clips, coins, keys, screws ... that could connect a terminal to another.
- Do not short-circuit the battery terminals as this may cause an explosion or fire.

- The battery/ies should be stored at a temperature of below 40°C.
- CAUTION: Do not use the appliance if the glass is cracked or broken.
- Use the device only with the battery/ies for which has been specifically designed. The use of any other battery/ies can cause an explosion or fire hazard.

ASSEMBLING THE BATTERY/S

- Caution: During the handling of batteries, do not touch both poles at the same time, as this will provoke the partial discharge of the stocked energy therefore affecting longevity.
- Remove the cover of the battery compartment.
- Check that the plastic covering that protects the battery has been removed (some batteries are sold with a protective covering).
- Put the battery in its compartment, respecting the polarity (Fig.1).
- Replace the cover of the battery compartment.
- It is essential that the batteries should always be the same kind and the same charge. Never mix alkaline batteries with normal one (Carbon-Zinc) and/or rechargeable ones.
- N.B. When changing the batteries, the configuration of the appliance will be removed, and you will need to reset the parameters of the appliance.

INSTRUCTIONS FOR USE

BEFORE USE:

- Make sure that all products' packaging has been removed.

WEIGHING UNITS SELECTOR KG/LB/ST:

- Set the unit of measure in "LB" or "KG" by pressing the button on the bottom side of your scale while the power is off. (Fig.2).
- After selecting the preferred measurement unit the selected measurement unit will be activated when the scale is turned back on. (Fig.3).

QUICKLY WEIGHING YOURSELF:

- If you would like to weigh yourself without going through the more advanced settings options, follow the below steps:
- Make sure the battery is installed and the scale is placed on a hard, even surface.

- Turn the appliance on, by pressing lightly on the weighing platform and wait for the screen to show 0.0.
- Place yourself in the centre of the platform avoiding abrupt movements, maintaining immobile and without touching other objects, until the appliance shows your weight.
- To avoid inaccurate results caused by electromagnetic interference between electrical and electronic equipment's, do not use this appliance near to a cell phone or microwave oven.
- The LCD will begin to display your weight. When the digits flash, the weight is stabilized and locked in.

SETTING UP YOUR PERSONAL DATA:

- To set up your personal data, follow these steps:
 1. With the scale o, press the "SET" button to power it on.
 2. With the scale on, you can now select a user profile. The scale's memory supports up to 13 users.
 3. To select a user profile, press the "SET" button and use the up and down arrow buttons to navigate through each user profile from "p-01" to "p-13."
 4. With a profile selected, you can now set your gender. Use the up and down arrow buttons located next to the "SET" button to switch between male or female, indicated by a man or woman figure on the left side of the screen. Press "SET" again to lock in the gender and move to the next setting.
 5. You can now set up your age. Use the up and down arrow buttons to adjust the number on screen to match your age. When you're done, press the "SET" button again to confirm your age and move on to the next setting.
 6. Finally, you can enter your height. Use the up and down arrow buttons to adjust the number on screen to match your height. Press "SET" again once you've confirmed your height to lock it in.
 7. The screen should now be showing 0.0, and is ready for you to weigh yourself.

MEASURING WITH THE SCALE:

- To obtain your measurements with the scale, follow these steps:
 1. With the scale off, press the "SET" button to power it on.
 2. Press "SET" again to enter user profile selection mode, and select your desired profile.
 3. When the screen on the scale displays "0.0", step onto the scale gently.
 4. If you would like your full set of measurements, you should step onto the scale bare feet, with nothing covering your feet, and with your feet making contact with each of the four metal contact points on the top of the scale.
 5. The scale will first measure your weight.
 6. The weight reading on the screen will stabilize and lock in, flashing once to indicate that the reading is locked.
 7. Next, the scale will gather the rest of your measurements and display them, showing body fat percentage, hydration percentage, muscle percentage, bone mass percentage, calories and BMI (Body Mass Index).
 8. The scale will display all of your measurements again before the screen turns off.

BODY FAT:

	Age	Underfat	Healthy	Overfat	Obese
	<40	<8%	8%-21%	22%-25%	>25%
	40-60	<11%	11%-22%	23%-28%	>28%
	>60	<13%	13%-25%	26%-30%	>30%
	<40	<20%	20%-33%	34%-39%	>39%
	40-60	<21%	21%-34%	35%-40%	>40%
	>60	<22%	22%-36%	37%-42%	>42%

BODY WATER

	Age	Underfat	Healthy	Overfat
	<40	<50%	50%-65%	>65%
	40-60	<48.5%	48.5%-63.5%	>63.5%
	>60	<47%	47%-62%	>62%
	<40	<45%	45-60%	>60%
	40-60	<43.5%	43.5%-58.5%	>58.5%
	>60	<42%	42-57%	>57%

MUSCLE MASS

	Age	Underfat	Healthy	Overfat
	<40	<40%	40%-53%	>53%
	40-60	<38%	38%-52%	>52%
	>60	<37%	37%-50%	>50%
	<40	<32%	32%-43%	>43%
	40-60	<28%	28%-41%	>41%
	>60	<28%	28%-40%	>40%

BONES MASS:

	Underfat	Healthy	Overfat
	<12%	12%-18%	>18%
	<10%	10%-16%	>16%

BASAL METABOLIC INDEX:

Indicates the number of calories the body needs to perform basic metabolic functions.

Underfat	Healthy
<1500	>1500

BMI (BODY MASS INDEX):

BMI=kg/m²

Low	Healthy	High	Very High
Underfat	Healthy	Overfat	Obese
<18.5	18.5-24.9	25.0-29.9	>30

In the event of doubt, consult your doctor or pharmacist for clarification.

STAND-BY FUNCTION:

- In order to save energy, the appliance passes to stand-by after few seconds if the user does not use it.
- To return to normal operation simply press lightly on the weighing platform and wait for the screen to show 0,0.

ONCE YOU HAVE FINISHED USING THE APPLIANCE:

- The device turns off automatically after the use.
- Remove the batteries from the appliance if you are not intending to use the appliance during an extended period.
- Clean the appliance.

SPECIAL MESSAGES:

- The following messages can appear on the screen:
- Err: means that the appliance's capacity has been exceeded or that there is an error and it is necessary to re-test.
- ErrL: means that there is an error in the body fat percentage measurement. It is very low.
- ErrH: means that there is an error in the body fat percentage measurement. It is very high.
- Lo: means that the battery has to be changed.

CLEANING

- Clean the equipment with a damp cloth with a few drops of washing-up liquid and then dry.
- Do not use solvents, or products with an acid or base pH such as bleach, or abrasive products, for cleaning the appliance.
- Do not let water or any other liquid get into the air vents to avoid damage to the inner parts of the appliance.
- Never submerge the appliance in water or any other liquid or place it under running water.

SUPPLIES

- Always use original supplies, designed specifically for your appliance model.
- Only the following type of supplies must be used with the appliance. (Fig.4).
- Batteries: 2 AAA.
- These supplies are available in specialist shops.

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English

WARRANTY AND TECHNICAL ASSISTANCE

This product enjoys the recognition and protection of the legal guarantee in accordance with current legislation. To enforce your rights or interests you must go to any of our official technical assistance services.

You can find the closest one by accessing the following web link: <http://taurus-home.com/>

You can also request related information by contacting us.

You can download this instruction manual and its updates at <http://taurus-home.com/>

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